

AirTec

alternating air system

The AirTec alternating air system mattresses provide advanced clinical benefits in pressure redistribution and in the treatment of pressure sores, while offering excellent user comfort and support.



AirTec 5 - Alternating Pressure Mattress Overlay

Order code: M20-AMP (S/KS)

The AirTec 5 – 5” overlay provides cost effective care for users with low to medium risk of pressure sores, up to stage II.

Dimensions: Single 1980 x 880 x 130mm, King Single 1980 x 1050 x 130mm

Features

- Single and King Single
- Whisper quiet digital pump
- A & B cells alternating mode - continuous operation
- Static/nursing operation mode
- Alternate or Static seat inflate mode
- Adjustable cycle times
- Adjustable support weight and comfort settings
- Low pressure alarm with audible buzzer and indicator
- Power failure and service support alarm
- Smart pressure detection technology - prevents the mattress from being over pressurised
- Quick connect plugs
- Full zip-around 2-way stretch Nylon TPU cover

AirTec 7.5 - Alternating Pressure Mattress

Order code: M21-AMP (S/KS)

The AirTec 7.5 – 7.5” full mattress replacement features 5” alternating air cells with 2.5” foam underlay. Suitable for users with medium risk/stage III pressure sores.

Dimensions: Single 1980 x 880 x 190mm, King Single 1980 x 1050 x 190mm



cobalthealth.com.au

sales@cobalthealth.com.au
+61 2 8853 1100

Unit 34, 9 Salisbury Rd, Castle Hill
NSW 2154 Australia